

# ML in Medical Imaging. Multi-Task Learning

---

## 1. Downloading the data

---

1. Clone/download the code file
2. Download the image dataset from <https://zenodo.org/record/7013610#.ZEgRR0zP23I>. Store the images in a folder named `data` **at the same level** as the main folder that contains the code.
3. (Optional) Download the trained models from [https://liveuclac-my.sharepoint.com/:f:/g/personal/rmapasr\\_ucl\\_ac\\_uk/EldJqhPlzvhDjnJZb0GiYYABB\\_NJhURn98r4nxJ55A7sxw?e=M7L7fg](https://liveuclac-my.sharepoint.com/:f:/g/personal/rmapasr_ucl_ac_uk/EldJqhPlzvhDjnJZb0GiYYABB_NJhURn98r4nxJ55A7sxw?e=M7L7fg) Store them in a folder named `models` **inside** the main folder that contains the code.

## 2. Set up the environment

---

Run the following script to set up a new Conda environment and install the necessary packages:

```
conda create --name mphy0041-cw2 -c conda-forge numpy nibabel matplotlib tensorflow=2.10 pytorch=1.12 monai
conda install -c anaconda pathlib
```

## 3. Train & test the models

---

Modify the first few lines of code in the file `main.py`. Each parameter can take a value of 1 (True) or 0 (False):

```
# Choose whether to train and/or test model(s)
TRAIN          = 1
TEST           = 1

# Choose which models to test
BASE_CASE      = 1
AUX_SEGMENT_3  = 1
AUX_SEGMENT_6  = 1
AUX_RECONSTRUCT = 1
```

Run the following script on the command line:

```
python main.py
```